

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Ethylamine - 99.0%

Product Number: 724033

CAS Number: 75-04-7

Formula: C₂H₇N

Formula Weight: 45.08 g/mol

TEST	Specification
Appearance (Color)	Colorless to Blue
Colorless to Blue (if liquid phase)	
Colorless (if gas phase)	
Appearance (Form)	Gas or Liquid
Infrared spectrum	Conforms to Structure
GC (area %)	≥ 99.0 %

Specification: PRD.0.ZQ5.10000099454

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.