

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
5,7-Dichloro-1-indanone - 97%

Product Number: **692107**
CAS Number: 448193-94-0

Formula: C₉H₆Cl₂O
Formula Weight: 201.05 g/mol



TEST

Specification

Appearance (Color)
Appearance (Form)
Solid, Powder, Crystals, Crystalline Powder and/or Chunks
Infrared spectrum
Titration by AgNO₃
Purity (GC)

White to Beige and Faint Brown to Brown
Conforms to Requirements

Conforms to Structure
96.5 - 103.5 %
≥ 96.5 %

Specification: PRD.0.ZQ5.10000012610

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.