

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Silver chloride - ReagentPlus®, 99%

Product Number: 227927
CAS Number: 7783-90-6
MDL: MFCD00003399
Formula: AgCl
Formula Weight: 143.32 g/mol

AgCl

TEST**Specification**

Appearance (Color)
Appearance (Form)
X-Ray Diffraction
Titration by AgNO₃
 % Ag After NaCN

White to Off-White
Powder or Chunks
Conforms to Structure
74.1 - 76.3 %

Specification: PRD.0.ZQ5.10000026714

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.