

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Anti-Transglutaminase Type II antibody produced in goat - whole antiserum

Product Number: **T7066**

MDL: MFCD01324001

Storage Temperature: -20 °C

TEST	Specification
------	---------------

Appearance	Conforms
antiserum containing 0.05% sodium azide	
working dilution by immunoblot	Conforms
1:1000 to 1:2000 using mouse liver tissue extracts, anti-rabbit IGG peroxidase conjugate and Chemiluminescent detection system	
working dilution by immunoprecipitation	Conforms
10 uL is used to precipitate type II transglutaminase from Mouse liver and spleen tissue extracts	

Specification Date : 11/18/2010

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.