


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Tetrabutylammonium hexachloroplatinate(IV)

Product Number: 342149
CAS Number: 18129-78-7
MDL: MFCD00075248
Formula: C₃₂H₇₂Cl₆N₂Pt
Formula Weight: 892.72 g/mol



TEST	Specification
Appearance (Color)	Orange
Appearance (Form)	Powder or Crystals
Infrared spectrum	Conforms to Structure
Carbon	41.9 - 44.1 %
Nitrogen	2.9 - 3.3 %
ICP: Confirms Platinum Compone	Confirmed

Specification: PRD.0.ZQ5.10000027436

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.