


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
4-Oxo-TEMPO-d₁₆, ¹⁵N, free radical - 98 atom % ¹⁵N, 98 atom % D

Product Number: **487740**
CAS Number: 80404-14-4
MDL: MFCD00190575
Formula: C₉D₁₆¹⁵N₁O₂
Formula Weight: 186.98 g/mol



TEST	Specification
Appearance (Color)	Orange
Appearance (Form)	Solid
Purity (GC)	≥ 99.0 %
¹⁵ N Enrichment	≥ 98.0
Isotopic Enrichment-D	≥ 98.0
Melting Point	Conforms

Specification: PRD.0.ZQ5.10000103249

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

