

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Lead - shot, 1-3 mm, 99.995% trace metals basis

Product Number:

695912

Pb

CAS Number:

7439-92-1

Formula:

Pb

Formula Weight:

207.2 g/mol

TEST	Specification
Appearance (Color)	Conforms to Requirements
Grey	
Appearance (Form)	Conforms
Metal Shot	
Size	Conforms
1-3 mm	
ICP Major Analysis	Confirmed
Confirms Lead Component	
Trace Metal Analysis	≤ 55.0 ppm
Purity	Meets Requirements
99.995% Based On Trace Metals Analysis	

Specification: PRD.0.ZQ5.10000035434

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.