

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com**Product Specification**

Product Name:

Chemiluminescent Peroxidase Substrate-3 sufficient for 500 mL substrate

Product Number:

CPS3500

TEST	Specification
Appearance (Turbidity)	Clear
Appearance (Color)	Colorless
Appearance (Form)	Solution
Dot Blot Sensitivity	Conforms
Successfully detects a minimum of 5 ng of mouse IgG. Normal saline alone gives no significant signal and there is no excessive non-specific binding.	

Remarks:

Specification Date : 10/08/2010

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.