

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
(S,R,S)-AHPC-PEG₆-Azide - ≥95%Product Number: **909378**
CAS Number: 2086298-71-5Formula: C₃₆H₅₅N₇O₁₀S
Formula Weight: 777.94 g/mol
Storage Temperature: 2 - 8 °C

TEST	Specification
------	---------------

Appearance (Color)	White to Yellow to Brown
Appearance (Form)	Powder or Crystal or Chunk(s) or Solid or Liquid
May be gummy	
1H NMR Spectrum	Conforms to Structure
Carbon Content	52.8 - 58.4 %
Nitrogen Content	11.9 - 13.2 %
Purity (HPLC)	≥ 95.0 %
Specific Rotation	
Specification: Report Only	

Specification: PRD.0.ZQ5.10000101312

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.