

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Titanium(IV) oxide - ReagentPlus®, ≥99%

Product Number: 14021
CAS Number: 13463-67-7
MDL: MFCD00011269
Formula: TiO₂
Formula Weight: 79.87 g/mol



TEST	Specification
Arsenic (As)	≤ 0.002 %
Iron (Fe)	≤ 0.005 %
Lead (Pb)	≤ 0.005 %
Titanium Dioxide	≥ 99.0 %
Loss on Ignition (1000°C, 2 h)	≤ 0.5 %

Specification: PRD.0.ZQ5.10000103174

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

