

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Palladium(II) sulfate - 98%

Product Number: 205788
CAS Number: 13566-03-5
MDL: MFCD00011173
Formula: O4PdS
Formula Weight: 202.48 g/mol



TEST	Specification
------	---------------

Appearance (Color)	Red to Very Dark Red and Red-Brown and Very Dark Red-Brown
Appearance (Form)	Powder
ICP Major Analysis Confirms Pd and S Components	Confirmed
Palladium (Pd)	51.2 - 53.9 %

Specification: PRD.0.ZQ5.10000035693

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.