

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com


Product Specification

Product Name:

Hydroxylamine sulfate - 99.999% trace metals basis

Product Number:

379913



CAS Number:

10039-54-0

MDL:

MFCD00044869

Formula:

H6N2O2·H2O4S

Formula Weight:

164.14 g/mol

TEST	Specification
------	---------------

Appearance (Color)

White

Appearance (Form)

Conforms to Requirements

Crystals or Crystalline Powder

Infrared spectrum

Conforms to Structure

Titration by KMNO4

39.2 - 41.7 %

% H2NOH

ICP Major Analysis

Confirmed

Confirms Sulfur Component

Trace Metal Analysis

≤ 20.0 ppm

Purity

Conforms

99.999% Based On Trace Metals Analysis

Specification: PRD.0.ZQ5.10000030718

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.