

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Anti-Human IgG (γ -chain specific) antibody produced in goat - whole antiserum

Product Number:

I1136

Storage Temperature:

-20 °C

TEST**Specification**

Appearance (Form)	Liquid
Appearance (Turbidity)	Clear to Slightly Hazy
Appearance (Color)	Brown-Yellow
Protein by Biuret	40 - 100 mg/ml
Antibody Concentration(Method)	Quantitative Precipitin Assay
Antibody Concentration	4.0 - 6.0 mg/ml
ELISA - Titer	1:40,000-1:80,000
Recommended Retest Period	-----
5 Years	

Specification: PRE.0.ZQ5.10000019548

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.