


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Vanadium(III) oxide - 98%

Product Number: 215988
CAS Number: 1314-34-7
MDL: MFCD00011455
Formula: O₃V₂
Formula Weight: 149.88 g/mol

**TEST****Specification**

Appearance (Color)	Black
Appearance (Form)	Powder or Powder with Chunk(s)
Titration by KMNO ₄ (% V)	66.3 - 69.7 %
X-Ray Diffraction	Conforms to Structure

Specification: PRD.0.ZQ5.10000092048

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.