

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Silver chloride - anhydrous, beads, -10 mesh, 99.998% trace metals basis

Product Number:**449571**

CAS Number:

7783-90-6

MDL:

MFCD00003399

Formula:

AgCl

Formula Weight:

143.32 g/mol

AgCl

TEST	Specification
------	---------------

Appearance (Color)

Conforms to Requirements

White to Grey or Tan

Appearance (Form)

Beads

ICP Major Analysis

Confirmed

Confirms Silver Component

Purity

Conforms

99.998% Purity Based On Trace Metals Analysis

Trace Metal Analysis

≤ 25.0 ppm

Specification: PRD.2.ZQ5.10000027697

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.