

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Bongkreikic acid solution - from *Pseudomonas cocovenenans*, ≥95% (HPLC), ~1 mg/mL

Product Number:

B6179



Formula:

C₂₈H₃₈O₇

Formula Weight:

486.60 g/mol

Storage Temperature:

-20 °C

TEST	Specification
------	---------------

Appearance (Form)

Clear Solution

Conforms

Appearance (Colour)

Colorless to Light Yellow

Appearance (Solvent)

0.01M TRIS

Purity (HPLC)

≥ 95.00 %

Vial Content (Size)

500.00 - 500.00

Vial Content (Unit Of Measure)

µg

Vial Content (Overpack)

≥ 1.00 %

Recommended Retest Period

2 Years

Specification: PRE.1.ZQ5.10000001279

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.