


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
2,5-Dichloro-3-nitrobenzoic acid - 98%

Product Number: 429848
CAS Number: 88-86-8
MDL: MFCD00014691
Formula: C₇H₃Cl₂NO₄
Formula Weight: 236.01 g/mol



TEST

Specification

Appearance (Color)

Conforms to Requirements

White to Yellow to Tan

Appearance (Form)

Conforms to Requirements

Powder, Crystals, Crystalline Powder and/or Chunks

Infrared spectrum

Conforms to Structure

Purity (Titration by NaOH)

97.5 - 102.5 %

GC (area %)

≥ 97.5 %

Specification: PRD.0.ZQ5.10000036366

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.