


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Ruthenium(IV) oxide hydrate - 99.9% trace metals basis

Product Number: 463760
CAS Number: 32740-79-7
MDL: MFCD00149846
Formula: O₂Ru · xH₂O
Formula Weight: 133.07 g/mol



TEST	Specification
Appearance (Color)	Black
Appearance (Form)	Powder
ICP Major Analysis Confirms Ruthenium Component	Confirmed
Purity 99.9% Based On Trace Metals Analysis	Meets Requirements
Trace Metal Analysis	≤ 1500.0 ppm

Specification: PRD.0.ZQ5.10000072180

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.