

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com


Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

Product Number: F8257
CAS Number: 30484-77-6
MDL: MFCD00058198
Formula: C₂₆H₂₆F₂N₂ · 2HCl
Formula Weight: 477.42 g/mol
Storage Temperature: 2 - 8 °C



TEST

Specification

Appearance (Color)	White
Appearance (Form)	Powder
Infrared spectrum	Conforms to Structure
Solubility (Color)	Colorless
Solubility (Turbidity)	Clear
50 mg/mL, CHCl ₃ :MeOH (1:1)	
Water (by Karl Fischer)	≤ 2.0 %
Purity (TLC)	≥ 99 %

Specification: PRD.0.ZQ5.10000043946

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.