


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
3-Indoleglyoxylic acid - 98%

Product Number: 220019
CAS Number: 1477-49-2
MDL: MFCD00005625
Formula: C₁₀H₇NO₃
Formula Weight: 189.17 g/mol
Storage Temperature: 2 - 8 °C



TEST	Specification
Appearance (Color)	Conforms to Requirements
Yellow to Yellow-Orange to Gold	
Appearance (Form)	Conforms to Requirements
Powder, Crystals, Crystalline Powder and/or Chunks	
Infrared spectrum	Conforms to Structure
Purity (Titration by NaOH)	97.5 - 102.5 %
Purity (HPLC)	≥ 97.5 %

Specification: PRD.1.ZQ5.10000020377

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.