

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com


Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
(R)-(-)-2-Amino-1-pentanol - 97%

Product Number: 534579
CAS Number: 80696-30-6
MDL: MFCD02683228
Formula: C₅H₁₃NO
Formula Weight: 103.16 g/mol



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Conforms to Requirements
Solid, Powder or Crystalline Solid	
Infrared Spectrum	Conforms to Structure
Optical Rotation	
C = 1%, CHCl ₃	
Purity (GC)	≥ 97.5 %

Specification: PRD.0.ZQ5.10000100562

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.