

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com


Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
cis-3-Hexenyl butyrate - ≥98%, FG

Product Number: W340200
CAS Number: 16491-36-4
MDL: MFCD00036540
Formula: C10H18O2
Formula Weight: 170.25 g/mol



| TEST | Specification |
|----------------------------|-----------------------|
| Appearance (Color) | Colorless |
| Appearance (Form) | Liquid |
| Refractive index at 20 ° C | 1.429 - 1.433 |
| Infrared spectrum | Conforms to Structure |
| Purity (GC) | ≥ 98.0 % |
| Arsenic (As) | ≤ 3 ppm |
| Cadmium (Cd) | ≤ 1 ppm |
| Mercury (Hg) | ≤ 1 ppm |
| Lead (Pb) | ≤ 10 ppm |
| Expiration Date Period | ----- |
| 5 Years | |

Specification: PRD.4.ZQ5.1000004826

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.