

sigma-aldrich.com

3050 Spruce Street, Saint Louis, MO 63103, USA Website: www.sigmaaldrich.com Email USA: techserv@sial.com Outside USA: eurtechserv@sial.com

## **Product Specification**

Product Name: Scandium(III) oxide - 99.9% trace rare earth metals basis

Product Number:	307874
CAS Number:	12060-08-1
MDL:	MFCD00011223
Formula:	O3Sc2
Formula Weight:	137.91 g/mol



TEST	Specification	
Appearance (Color)	White	
Appearance (Form)	Powder	
Complexiometric EDTA	62.0 - 68.0 %	
% Scandium		
ICP Major Analysis	Confirmed	
Confirms Scandium Component		
Purity	Meets Requirements	
99.9% Based on Trace Rare Earth Metals Analysis		
Trace Rare Earth Analysis	<pre>&lt; 1500.0 ppm</pre>	

Specification: PRD.3.ZQ5.10000045940

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.