

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

LITHIUM, INGOT, 5.7CM DIAM., 99.9%

Product Number:**265977**

CAS Number:

7439-93-2

MDL:

MFCD00134051

Formula:

Li

Formula Weight:

6.94 g/mol

TEST	Specification
------	---------------

Appearance (Color)

Conforms to Requirements

Grey

Appearance (Form)

Conforms

Metal: Approximately 5.7cm Diameter

ICP Major Analysis

Confirmed

Confirms Lithium Component

Purity

Meets Requirements

99.9% Based On Trace Metals Analysis of ICP and AA combined

Trace Metal Analysis

≤ 1500.0 ppm

Atomic Absorption

≤ 1500.0 ppm

Specification: PRD.0.ZQ5.10000039596

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.