


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Hexanoic anhydride - 97%

Product Number: 194530
CAS Number: 2051-49-2
MDL: MFCD00009509
Formula: C₁₂H₂₂O₃
Formula Weight: 214.3 g/mol



TEST	Specification
------	---------------

Appearance (Clarity)	Clear
Appearance (Color)	Colorless to Very Faint Yellow
Appearance (Form)	Liquid
Refractive index at 20 ° C	1.428 - 1.430
Infrared spectrum	Conforms to Structure
Morpholine Titration by HCl	97.0 - 102.0 %

Specification: PRD.1.ZQ5.10000029697

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.