


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
m-Toluidine - 99%

Product Number: 511218
CAS Number: 108-44-1
MDL: MFCD00007808
Formula: C₇H₉N
Formula Weight: 107.15 g/mol



TEST

Specification

Appearance (Color)
Appearance (Form)
Infrared spectrum
Purity (GC)

Very Faint Yellow to Yellow and Very Faint Brown to Brown
Liquid
Conforms to Structure
≥ 98.5 %

Specification: PRD.0.ZQ5.10000028788

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.