


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Acridine - 97%

Product Number: A23609
CAS Number: 260-94-6
MDL: MFCD00005025
Formula: C₁₃H₉N
Formula Weight: 179.22 g/mol



TEST

Specification

Appearance (Color)

Light Yellow to Yellow and Light Beige

Appearance (Form)

to Beige and Light Brown to Brown and Yellow-Beige and Yellow-Brown

Infrared spectrum

Powder or Crystals

Purity (Titration by HCL04)

Conforms to Structure

Purity (HPLC)

96.5 - 103.5 %

≥ 96.5 %

Remarks:

Reviewing Chemist :

Specification: PRD.0.ZQ5.10000012163

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.