

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Sodium pyridine-2-sulfinate - ≥95%

Product Number: **900757**
CAS Number: 24367-66-6

Formula: C₅H₄NNaO₂S
Formula Weight: 165,14 g/mol



TEST	Specification
Appearance (Color)	White to Yellow
Appearance (Form)	Powder or Crystals
¹ H NMR Spectrum	Conforms to Structure
Carbon Content	34.5 - 38.2 %
Nitrogen Content	8.0 - 9.0 %
Sulfur Content	18.4 - 20.4 %
ICP Major Analysis	Confirmed
Confirms Sodium Component	

Specification: PRD.0.ZQ5.10000095135

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.