

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Strontium titanate – single crystal substrate,

Product Number: 634689
CAS Number: 12060-59-2
MDL: MFCD00049554
Formula: O3SrTi
Formula Weight: 183.49 g/mol



TEST	Specification
Appearance (Color)	White to Off-White
Appearance (Form)	Crystals
Particle Size	Conforms
10 mm x 10 mm x 0.5 mm	
ICP Major Analysis	Confirmed
Confirms Sr & Ti Components	
Orientation 100	Confirmed

Specification: PRD.0.ZQ5.10000040692

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.