

3050 Spruce Street, Saint Louis, MO 63103, USA
Website: www.sigmaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Antifoam 204

Product Number: **A6426**

TEST	Specification
Appearance (Form)	Viscous Liquid
Suitable For Use As A Defoamer	Suitable
Cloud Point 1% Solution	18.0- 21.0 °C
pH c = 2.5% in MeOH-H ₂ O (80-20)	5.00- 7.50
Color Test	≤70 APHA
Specific Gravity	1.01- 1.03
Recommended Retest Period 2 Years	-----

Remarks:

Specification Date : 09/10/2010

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.