

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
E-Toxate™ Water – endotoxin, free

Product Number: 2107 H₂O
CAS Number: 7732-18-5
MDL: MFCD00011332
Formula: H₂O
Formula Weight: 18.02 g/mol

| TEST | Specification |
|---------------------------------|-------------------------|
| Appearance (Turbidity) | Clear |
| Appearance (Color) | Colorless |
| Appearance (Form) | Liquid |
| DNase, Exonuclease Detection: | None Detected |
| NICKase, Endonuclease Detection | None Detected |
| Endotoxin Level | Less than 0.1 EU per mL |

Specification: PRD.0.ZQ5.10000026321

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.