


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Bis(cyclopentadienyl)cobalt(II)

Product Number: 339164
CAS Number: 1277-43-6
MDL: MFCD00013749
Formula: C₁₀H₁₀Co
Formula Weight: 189.12 g/mol
Storage Temperature: 2 - 8 °C



TEST

Specification

Appearance (Color)	Black
Appearance (Form)	Conforms to Requirements
Crystals or Crystalline Powder and/or Chunks	
Infrared spectrum	Conforms to Structure
Carbon	60.0 - 67.0 %

Specification: PRD.0.ZQ5.10000021239

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.