

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.com


Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
2,3,7,8,12,13,17,18-Octaethyl-21H,23H-porphine zinc(II) - 97%

Product Number: 258423
CAS Number: 17632-18-7
MDL: MFCD00012156
Formula: C₃₆H₄₄N₄Zn
Formula Weight: 598.15 g/mol



TEST	Specification
Appearance (Color)	Conforms to Requirements
Dark Red to Purple	
Appearance (Form)	Conforms to Requirements
Powder, Fibers or Crystals	
Infrared spectrum	Conforms to Structure
Wavelength	399 - 405 nm
0.001-0.002g/l; Toluene	
Extinction Coefficient	≥ 200000
Purity (HPLC)	≥ 96.5 %

Specification: PRD.0.ZQ5.10000056089

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.