


3050 Spruce Street, Saint Louis, MO 63103, USA
Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
APhos Pd G4 - 95%

Product Number: 900608

Formula: C₃₀H₄₃N₂O₃PPdS
Formula Weight: 649.13 g/mol



TEST

Specification

Appearance (Color)	White to Beige
Appearance (Form)	Powder or Crystals
¹ H NMR Spectrum	Conforms to Structure
³¹ P NMR Spectrum	Conforms to Structure
Carbon Content	52.4 - 58.6 %
Nitrogen Content	4.1 - 4.6 %
Sulfur Content	4.6 - 5.2 %
ICP Major Analysis	Confirmed
Confirms Palladium Component	

Specification: PRD.1.ZQ5.10000102909

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.