

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
L-Aspartic acid-2-¹³C,¹⁵N - 98 atom % ¹⁵N, 99 atom % ¹³C

Product Number: **607703**

MDL: MFCD00083896



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder
Water (by Karl Fischer)	≤ 0.50 %
Purity (HPLC)	≥ 99.0 %
Chiral Purity	≥ 99 %
¹ H NMR Spectrum	Conforms to Structure
¹³ C NMR Spectrum	Conforms to Structure
¹³ C Enrichment	≥ 99.0 atom %
¹⁵ N Enrichment	≥ 98 atom %

Specification: PRD.0.ZQ5.10000102364

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.