


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
1,4-Butanediol dimethacrylate - 95%, contains 200-300 ppm MEHQ as inhibitor

Product Number: 234958
CAS Number: 2082-81-7
MDL: MFCD00008594
Formula: C₁₂H₁₈O₄
Formula Weight: 226.27 g/mol
Storage Temperature: 2 - 8 °C



TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Infrared spectrum	Conforms to Structure
Purity (GC)	≥ 94.5 %
MEHQ as Inhibitor	200 - 300 ppm

Specification: PRD.0.ZQ5.10000072211

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.