


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
(Acetylacetonato)dicarbonylrhodium(I) - 98%

Product Number: 288101
CAS Number: 14874-82-9
MDL: MFCD00009884
Formula: C₇H₇O₄Rh
Formula Weight: 258.03 g/mol



TEST	Specification
Appearance (Color)	Conforms to Requirements
Green to Dark Green to Brown	
Appearance (Form)	Conforms to Requirements
Crystals or Needles or Powder	
Infrared spectrum	Conforms to Structure
Carbon	31.8 - 33.4 %

Specification: PRD.0.ZQ5.10000025578

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.