

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Sulfosuccinic acid solution - 70 wt. % in H₂O

Product Number: **338567**

MDL: MFCD00075139
Formula: C₄H₆O₇S
Formula Weight: 198.15 AMU



TEST

Specification

Appearance (Color)
Appearance (Form)
Infrared spectrum
Purity (Titration by NaOH)

Colorless
Viscous Liquid
Conforms to Structure
68.0 - 72.0 %

Specification Date : 04/18/2011

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.