

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Hydrocortisone 21-acetate - ≥98% (HPLC)

Product Number:

H4126

CAS Number:

50-03-3

MDL:


MFCD00037714

Formula:

C₂₃H₃₂O₆

Formula Weight:

404.50 g/mol



TEST

Specification

Appearance (Color)

White to Off-White

Appearance (Form)

Powder

UV/Vis absorbance

242 - 243

lambda max 242 to 243 nm, in Methanol

Purity (HPLC)

≥ 98 %

Recommended Retest Period

5 Years

Specification: PRD.1.ZQ5.10000050758

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.