


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
PCy₃ Pd G4

Product Number: 900287

Formula: C₃₂H₄₈NO₃PPdS
Formula Weight: 664.19 g/mol

TEST	Specification
------	---------------

Appearance (Color)	White to Very Dark Yellow to Very Dark Brown and Faint Grey to Very Dark Grey
Appearance (Form)	Powder or Crystals
ICP Major Analysis	Confirmed
Confirms Palladium Component	
1H NMR Spectrum	Conforms to Structure
31P NMR Spectrum	Conforms to Structure
Carbon Content	54.7 - 61.1 %
Nitrogen Content	
Specification: Report Only	
Sulfur Content	
Specification: Report Only	

Specification: PRD.1.ZQ5.10000079435

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.