

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Hemoglobin from bovine blood - lyophilized powder

Product Number: H2500
CAS Number: 9008-02-0
MDL: MFCD00131282

Storage Temperature: 2 - 8 °C

TEST	Specification
------	---------------

Appearance (Color)	Brown-Red to Very Dark Brown-Red
Appearance (Form)	Powder
Solubility (Color)	Dark Brown-Red
Solubility (Form)	Solution
20 mg/mL H ₂ O	
Iron (Fe)	0.25 - 0.35 %
Water (by Karl Fischer)	≤ 8 %

Specification: PRD.1.ZQ5.10000065120

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.