

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Potassium isopropyltrifluoroborate – 97%

Product Number:

667153

MDL:


MFCD04112718

Formula:

C₃H₇BF₃K

Formula Weight:

149.99 g/mol



TEST

Specification

Appearance (Color)

White to Beige to Brown

Appearance (Form)

Powder or Crystals

Infrared spectrum

Conforms to Structure

Carbon Content

23.1 - 24.9 %

1H NMR Spectrum

Conforms to Structure

19F NMR Spectrum

Conforms to Structure

Specification: PRD.2.ZQ5.10000005470

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.