

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Corticosterone-9,11,12,12-d4 - 98 atom % D, 97% (CP)

Product Number: **802905**

Formula: C₂₁D₄H₂₆O₄
Formula Weight: 350 g/mol
Storage Temperature: -20 °C



TEST	Specification
------	---------------

Appearance (Color)	White
Appearance (Form)	Powder
Purity (HPLC)	≥ 97 %
Proton NMR Spectrum	Conforms
Enrichment (atom% D)	≥ 98 %

Specification: PRD.0.ZQ5.10000059852

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.