


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
12-Aminododecanoic acid - 95%

Product Number: 159247
CAS Number: 693-57-2
MDL: MFCD00008153
Formula: C₁₂H₂₅NO₂
Formula Weight: 215.33 g/mol



TEST	Specification
Appearance (Color)	White to Off White
Appearance (Form)	Powder or Crystals
Infrared spectrum	Conforms to Structure
Titration with HClO ₄	94.0 - 106.0 %
Purity (TLC)	≥ 95 %
Solubility (Turbidity)	Clear to Slightly Hazy
50mg/ml (5%), Acetic Acid	
Solubility (Color)	Colorless to Light Yellow

Specification: PRD.0.ZQ5.10000085194

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.