

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

Sodium chloride - BioReagent, suitable for cell culture, suitable for insect cell culture, suitable for plant cell

culture, ≥99%

Product Number: S5886 CAS Number: 7647-14-5

MDL: MFCD00003477

Formula: CINa

Formula Weight: 58.44 g/mol

NaCl

TEST Specification

Appearance (Color) White

Appearance (Form) Powder or Crystals

Titration with AgNO3 99.0 - 101.5 %

Cell Culture Test Pass
Insect Cell Test Pass
Plant Cell Culture Test Pass

Recommended Retest Period

6 Years

Specification: PRD.3.ZQ5.10000024677

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.