


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
2',7'-Dichlorofluorescein - ACS reagent

Product Number: 410217
CAS Number: 76-54-0
MDL: MFCD00005047
Formula: C₂₀H₁₀Cl₂O₅
Formula Weight: 401.20 g/mol



TEST	Specification
Appearance (Color)	Conforms to Requirements
Orange to Red-Brown	
Appearance (Form)	Powder
Infrared spectrum	Conforms to Structure
Clarity of Solution	Pass
Use Test	Pass
Meets ACS Requirements	Conforms

Specification: PRD.0.ZQ5.10000072842

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.