


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
1,3-Diiminoisoindoline - 97%

Product Number: 279528
CAS Number: 3468-11-9
MDL: MFCD00010062
Formula: C₈H₇N₃
Formula Weight: 145.16 g/mol



TEST

Specification

Appearance (Color)	Off-White to Yellow
Appearance (Form)	Powder to Crystalline Powder
Proton NMR spectrum	Conforms to Structure
Purity (HPLC)	≥ 96.5 %
Titration with HClO ₄	96.5 - 103.5 %

Specification: PRD.1.ZQ5.10000053732

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.