

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Calcium nitrate tetrahydrate – suitable for cell culture, suitable for plant cell culture

Product Number:

C2786

CAS Number:

13477-34-4

MDL:


MFCD00149604

Formula:

CaN2O6 · 4H2O

Formula Weight:

236.15 g/mol



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder or Crystals
Solubility (Turbidity)	Clear
Solubility (Color)	Colorless
c = 20g/100mL; Water	
Complexometric EDTA	≥ 99.0 %
Plant Cell Culture Test	Pass
Cell Culture Test	Pass
Recommended Retest Period	-----
3 years	

Specification: PRD.0.ZQ5.10000019330

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.