

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

4-Aminoresorcinol hydrochloride - 96%

Product Number:

163171

CAS Number:

34781-86-7

MDL:


MFCD00012963

Formula:

C<sub>6</sub>H<sub>7</sub>NO<sub>2</sub> · HCl

Formula Weight:

161.59 g/mol



### TEST

### Specification

Appearance (Color)

Light Blue to Dark Blue and Light Purple  
to Dark Purple and Light Red to Dark Red and Light Grey to Very  
Dark Grey

Appearance (Form)

Powder

Proton NMR spectra

Conforms to Structure

Titration by AgNO<sub>3</sub>

95.5 - 104.5 %

Specification: PRD.2.ZQ5.10000033834

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.