

3050 Spruce Street, Saint Louis, MO 63103, USA

KH₂PO₄

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

Potassium phosphate monobasic - powder, suitable for cell culture, suitable for insect cell culture, suitable for plant cell culture, ≥99.0%

Product Number: P5655 CAS Number: 7778-77-0

MDL: MFCD00011401

Formula: H2KO4P Formula Weight: 136.09 g/mol

TEST Specification

Appearance (Color) White
Appearance (Form) Powder
Solubility (Color) Colorless
Solubility (Turbidity) Clear

100 mg/ml, H2O

Potassium Positive

ID by ICP

Phosphorus (P) Positive

ID by ICP

Titration with NaOH

Cell Culture Test

Pass
Insect Cell Test

Pass

Plant Cell Culture Test

Pass

Recommended Retest Period ------

4 years

Specification: PRD.1.ZQ5.10000016382

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.